

ARG80984
Dog CAV / Canine Adenovirus IgG antibody ELISA KitPackage: 96 wells
Store at: 4°C

Summary

Product Description	ARG80984 Dog CAV / Canine Adenovirus IgG antibody ELISA Kit is an Enzyme Immunoassay kit for the determination of Dog's CAV / Canine Adenovirus IgG antibody in serum and plasma.
Tested Reactivity	Dog
Tested Application	ELISA
Target Name	Dog CAV / Canine Adenovirus IgG
Conjugation	HRP
Conjugation Note	TMB at 450 nm
Sample Type	Serum and plasma.
Sample Volume	100 µl

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Microbiology and Infectious Disease kit
---------------	---