

ARG80985 Dog CCV / Canine Corona Virus IgG antibody ELISA Kit

Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG80985 Dog CCV / Canine Corona Virus IgG antibody ELISA Kit is an Enzyme Immunoassay kit for the determination of Dog's CCV / Canine Corona Virus IgG antibody in serum and plasma.
Tested Reactivity	Dog
Tested Application	ELISA
Target Name	Dog CCV / Canine Corona Virus IgG
Conjugation	HRP
Conjugation Note	TMB at 450 nm
Sample Type	Serum and plasma.
Sample Volume	100 µl

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Coronavirus antibodies ; Coronavirus ELISA Kits ; Coronavirus recombinant proteins ; Related news: ACE2, receptor of 2019-nCoV Exploring Antiviral Immune Response
Research Area	Microbiology and Infectious Disease kit