

ARG81004 Feline Corona Virus / FCoV antibody ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG81004 Feline Corona Virus / FCoV antibody ELISA Kit is an Enzyme Immunoassay kit for the determination of Cat's Corona Virus / FCoV antibody in serum, plasma and ascitic fluid.
Tested Reactivity	Cat
Tested Application	ELISA
Target Name	Corona Virus / FCoV antibody
Conjugation	HRP
Conjugation Note	TMB at 450 nm
Sample Type	Serum, plasma and ascitic fluid.
Sample Volume	100 µl

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Coronavirus antibodies ; Coronavirus ELISA Kits ; Coronavirus recombinant proteins ; Related news: ACE2, receptor of 2019-nCoV Exploring Antiviral Immune Response
Research Area	Microbiology and Infectious Disease kit