

ARG81017 Infectious Bovine Rhinotracheitis / BHV-1 antibody ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG81017 Infectious Bovine Rhinotracheitis / BHV-1 antibody ELISA Kit is an Enzyme Immunoassay kit for the quantification of Bovine's Infectious Bovine Rhinotracheitis / BHV-1 antibody in serum and plasma
Tested Reactivity	Bov
Tested Application	ELISA
Target Name	Infectious Bovine Rhinotracheitis / BHV-1
Conjugation	HRP
Conjugation Note	TMB at 450 nm
Sample Type	Serum and plasma
Sample Volume	50 µl

Application Instructions

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.