

## ARG81031 Coconut ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary


Product Description	ARG81031 Coconut ELISA Kit is an Enzyme Immunoassay kit for the quantification of Coconut in food (extraction, dilution, centrifugation).
Tested Reactivity	Other
Tested Application	ELISA
Target Name	Coconut
Conjugation	HRP
Conjugation Note	TMB at 450 nm
Sensitivity	0.4 ppm
Sample Type	Food (extraction, dilution, centrifugation).
Standard Range	2 - 30 ppm
Sample Volume	100 µl

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related products:  New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>
Research Area	Immune System kit



ARG81031 Coconut ELISA Kit standard curve image

ARG81031 Coconut ELISA Kit results of a typical standard run with optical density reading at 450 nm.