

## ARG81032 Macadamia nut ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG81032 Macadamia nut ELISA Kit is an Enzyme Immunoassay kit for the quantification of Macadamia nut in food (extraction, centrifugation, dilution).
Tested Reactivity	Other
Tested Application	ELISA
Target Name	Macadamia nut
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.1 ppm
Sample Type	Food (extraction, centrifugation, dilution).
Standard Range	1 - 40 ppm
Sample Volume	100 µl


### Application Instructions

#### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

#### Bioinformation

Resrarch Area	Immune System kit
---------------	-------------------



ARG81032 Macadamia nut ELISA Kit standard curve image

ARG81032 Macadamia nut ELISA Kit results of a typical standard run with optical density reading at 450 nm.

---