

## ARG81044 Human Pregnenolone ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary


Product Description	ARG81044 Human Pregnenolone ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Pregnenolone in serum.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Pregnenolone
Conjugation	HRP
Conjugation Note	TMB at 450 nm
Sensitivity	0.05 ng/ml
Sample Type	Serum
Standard Range	0.1 - 25.6 ng/ml
Sample Volume	50 µl

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Signaling Transduction kit
---------------	----------------------------



ARG81044 Human Pregnenolone ELISA Kit standard curve image

ARG81044 Human Pregnenolone ELISA Kit results of a typical standard run with optical density reading at 450 nm.