

ARG81173 Goat IgG ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

| | |
|---------------------|---|
| Product Description | ARG81173 Goat IgG ELISA Kit is an Enzyme Immunoassay kit for the quantification of Goat IgG in serum, plasma and cell culture supernatants. |
| Tested Reactivity | Goat |
| Tested Application | ELISA |
| Target Name | IgG |
| Conjugation | HRP |
| Conjugation Note | Substrate: TMB and read at 450 nm |
| Sensitivity | 3.9 ng/ml |
| Sample Type | Serum, plasma and cell culture supernatants. |
| Standard Range | 3.9 - 250 ng/ml |
| Sample Volume | 100 µl |
| Alternate Names | Immunoglobulin G |

Application Instructions


| | |
|------------|---------|
| Assay Time | 3 hours |
|------------|---------|

Properties

| | |
|---------------------|--|
| Form | 96 well |
| Storage instruction | Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |

Bioinformation

| | |
|-----------|--|
| Highlight | Related products: IgG antibodies ; IgG ELISA Kits ; IgG Duos / Panels ; IgG recombinant proteins ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData |
|-----------|--|



ARG81173 Goat IgG ELISA Kit standard curve image

ARG81173 Goat IgG ELISA Kit results of a typical standard run with optical density reading at 450 nm.