

ARG81226 BPDE Protein Adduct ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG81226 BPDE Protein Adduct ELISA Kit is an Enzyme Immunoassay kit for the quantification of BPDE Protein Adduct in purified protein, plasma, serum, cell lysate.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	BPDE Protein Adduct
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	10 µg/ml
Sample Type	Purified protein, plasma, serum, cell lysate
Standard Range	15.6 - 1000 ng/ml
Sample Volume	100 µl

Application Instructions


Assay Time	~ 6 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
-----------	----------------------------------------------------------------------------------------------------------------------



ARG81226 BPDE Protein Adduct ELISA Kit standard curve image

ARG81226 BPDE Protein Adduct ELISA Kit results of a typical standard run with optical density reading at 450 nm.