

## ARG81278 Mouse GM-CSF ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG81278 Mouse GM-CSF ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse GM-CSF in serum, plasma or cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Target Name	GM-CSF
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	16 pg/ml
Sample Type	Serum, plasma or cell culture supernatants.
Standard Range	31.25 - 2000 pg/ml
Sample Volume	100 µl
Alternate Names	GM-CSF; Granulocyte-macrophage colony-stimulating factor; CSF; GMCSF; Molgramostin; Colony-stimulating factor; Sargramostim

### Application Instructions

Assay Time	~ 3.5 hours
------------	-------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.


### Bioinformation

Gene Symbol	Csf2
Gene Full Name	colony stimulating factor 2 (granulocyte-macrophage)
Background	The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of granulocytes and macrophages. The active form of the protein is found extracellularly as a homodimer. This gene has been localized to a cluster of related genes at chromosome region 5q31, which is known to be associated with interstitial deletions in the 5q- syndrome and acute myelogenous leukemia. Other genes in the cluster include those encoding interleukins 4, 5, and 13. [provided by RefSeq, Jul 2008]
Function	Cytokine that stimulates the growth and differentiation of hematopoietic precursor cells from various lineages, including granulocytes, macrophages, eosinophils and erythrocytes. [UniProt]

Highlight

Related products:  
[GM-CSF antibodies](#); [GM-CSF ELISA Kits](#); [GM-CSF Duos / Panels](#); [GM-CSF recombinant proteins](#);  
Related news:  
[HMGB1 in inflammation](#)  
[Inflammatory Cytokines](#)  
New ELISA data calculation tool:  
[Simplify the ELISA analysis by GainData](#)

Images



ARG81278 Mouse GM-CSF ELISA Kit standard curve image

ARG81278 Mouse GM-CSF ELISA Kit results of a typical standard run with optical density reading at 450 nm.