

ARG82005 Free Fatty Acid Assay Kit

Package: 100 tests
Store at: -20°C

Summary

Product Description	ARG82005 Free Fatty Acid Assay Kit is a detection kit for the quantification of Free Fatty Acid in serum, plasma, urine, saliva, milk, cell cultures and food.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Free Fatty Acid
Sample Type	Serum, plasma, urine, saliva, milk, cell cultures and food.
Standard Range	Colorimetric Assay: 7 - 1000 μ M (read at 570 nm) Fluorimetric Assay: 7 - 100 μ M (read at Ex / Em: 530 / 585 nm)
Sample Volume	10 μ l


Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	30 min

Properties


Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images



ARG82005 Free Fatty Acid Assay Kit standard curve image

Colorimetric Assay: ARG82005 Free Fatty Acid Assay Kit results of a typical standard run with optical density reading at 570 nm.



ARG82005 Free Fatty Acid Assay Kit standard curve image

Fluorimetric Assay: ARG82005 Free Fatty Acid Assay Kit results of a typical standard run.