

ARG82026 Pectate Lyase Assay Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82026 Pectate Lyase Assay Kit is a detection kit for the quantification of Pectate Lyase Activity in tissue extracts and cell lysate.
Tested Reactivity	All
Tested Application	FuncSt
Target Name	Pectate Lyase
Conjugation Note	Read at 235 nm.
Sample Type	Tissue extracts and cell lysate.
Sample Volume	10 µl

Application Instructions

Assay Time	15 min
------------	--------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.