

ARG82030 Proline Assay Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82030 Proline Assay Kit is a detection kit for the quantification of Proline Content in urine, serum, plasma, tissue extracts, cell lysate and cell culture supernatants.
Tested Reactivity	All
Tested Application	FuncSt
Target Name	Proline
Conjugation Note	Read at 520 nm.
Detection Range	2 - 200 µg/ml
Sample Type	Urine, serum, plasma, tissue extracts, cell lysate and cell culture supernatants.
Sample Volume	50 µl


Application Instructions

Assay Time	30 min
------------	--------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images



ARG82030 Proline Assay Kit typical data demonstration image

ARG82030 Proline Assay Kit results of a typical data with optical density reading at 520 nm.