

ARG82032
Cysteine Assay KitPackage: 96 wells
Store at: 4°C

Summary


Product Description	ARG82032 Cysteine Assay Kit is a detection kit for the quantification of Cysteine Content in tissue extracts, cell lysate and cell culture supernatants.
Tested Reactivity	All
Tested Application	FuncSt
Target Name	Cysteine
Conjugation Note	Read at 600 nm.
Sensitivity	20 mmol/L
Detection Range	20 - 2000 mmol/L
Sample Type	Tissue extracts, cell lysate and cell culture supernatants.
Sample Volume	20 µl

Application Instructions

Assay Time	15 min
------------	--------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82032 Cysteine Assay Kit typical data demonstration image

ARG82032 Cysteine Assay Kit results of a typical data with optical density reading at 600 nm.