

**ARG82034**  
**Heme Assay Kit**Package: 96 wells  
Store at: 4°C

### Summary


Product Description	ARG82034 Heme Assay Kit is a detection kit for the quantification of Heme Content in blood, serum, plasma and urine.
Tested Reactivity	All
Tested Application	FuncSt
Target Name	Heme
Conjugation Note	Read at 505 nm.
Sensitivity	1 µmol/l
Detection Range	1 - 100 µmol/l
Sample Type	Blood, serum, plasma and urine.
Sample Volume	20 µl

### Application Instructions

Assay Time	10 min
------------	--------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82034 Heme Assay Kit typical data demonstration image

ARG82034 Heme Assay Kit results of a typical data with optical density reading at 505 nm.