

**ARG82041**  
**Pectinase Assay Kit**Package: 96 wells  
Store at: 4°C

### Summary


Product Description	ARG82041 Pectinase Assay Kit is a detection kit for the quantification of Pectinase Activity in tissue extracts, cell lysate and cell culture supernatants.
Tested Reactivity	All
Tested Application	FuncSt
Target Name	Pectinase
Conjugation Note	Read at 540 nm.
Sensitivity	0.1 mg/ml
Detection Range	0.1 - 2 mg/ml
Sample Type	Tissue extracts, cell lysate and cell culture supernatants.
Sample Volume	20 µl

### Application Instructions

Assay Time	45 min
------------	--------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82041 Pectinase Assay Kit typical data demonstration image

ARG82041 Pectinase Assay Kit results of a typical data with optical density reading at 540 nm.