

ARG82121 Creatinine Assay Kit

Package: 500 tests
Store at: 4°C, -20°C

Component

Cat. No.	Component Name	Package
ARG82121-001	Reagent A	50 ml
ARG82121-002	Reagent B	50 ml
ARG82121-003	Creatinine Standard (50 mg/dl)	1 ml

Summary


Product Description	ARG82121 Creatinine Assay Kit is a detection kit for the quantification of Creatinine in serum, plasma and urine.
Tested Reactivity	All
Tested Application	FuncSt
Target Name	Creatinine
Conjugation Note	Read at 510 nm.
Sensitivity	0.1 mg/dl (8 µM)
Sample Type	Serum, plasma and urine.
Standard Range	5 - 50 mg/dl
Sample Volume	Serum and plasma: 30 µl Urine: 5 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	20 min

Properties

Form	Liquid
Storage instruction	Store components at 4°C or -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82121 Creatinine Assay Kit standard curve image

ARG82121 Creatinine Assay Kit results of a typical standard run with optical density reading at 510 nm.