

## ARG82126 Acetate Assay Kit (Colorimetric)

Package: 100 tests  
Store at: -20°C

### Summary


Product Description	ARG82126 Acetate Assay Kit (Colorimetric) is a detection kit for the quantification of Acetate in serum, plasma and food.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Acetate
Conjugation Note	Read at 570 nm (colorimetric), or 530/585 nm (fluorometric)
Sensitivity	0.13 mM
Detection Range	0.20 - 20 mM (colorimetric); 0.13 - 2 mM (fluorometric)
Sample Type	serum, plasma and food
Sample Volume	10 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	30 min

### Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82126 Acetate Assay Kit (Colorimetric) typical data demonstration image

ARG82126 Acetate Assay Kit (Colorimetric) results of a typical data with optical density reading at 570 nm.