

ARG82131
ADP Assay Kit (Fluorometric)Package: 100 tests
Store at: -20°C

Summary


Product Description	ARG82131 ADP Assay Kit (Fluorometric) is a detection kit for the quantification of ADP in cell culture supernatants.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	ADP
Conjugation Note	Read at 530/590 nm (fluorometric)
Sensitivity	0.1 µM
Detection Range	6 - 20 µM
Sample Type	cell culture supernatants
Sample Volume	40 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	30 min

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82131 ADP Assay Kit (Fluorometric) typical data demonstration image

ARG82131 ADP Assay Kit (Fluorometric) results of a typical data.