

ARG82132
ADP / ATP Ratio Assay Kit (Chemiluminescent)Package: 100 tests
Store at: -20°C

Summary

Product Description	ARG82132 ADP / ATP Ratio Assay Kit (Chemiluminescent) is a detection kit for the quantification of ADP / ATP Ratio in cells.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	ADP / ATP Ratio
Conjugation Note	Detect signals by Luminometer.
Sample Type	cells
Sample Volume	10 µl cultured cells (10 ³ - 10 ⁴)

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	20 min

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.