

ARG82153 Formaldehyde Assay Kit (Fluorometric)

Package: 100 tests
Store at: 4°C

Summary


Product Description	ARG82153 Formaldehyde Assay Kit (Fluorometric) is a detection kit for the quantification of Formaldehyde in biological, food, beverage and environmental samples.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Formaldehyde
Conjugation Note	Read at 370/470 nm (fluorometric)
Sensitivity	1.5 µM (45 ppb)
Detection Range	5 - 100 µM
Sample Type	biological, food, beverage and environmental samples
Sample Volume	50 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	30 min

Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82153 Formaldehyde Assay Kit (Fluorometric) typical data demonstration image

ARG82153 Formaldehyde Assay Kit (Fluorometric) results of a typical data.