

ARG82173
Ketone Body Assay Kit (Colorimetric)

Package: 200 tests

Store at: -20°C

Summary


Product Description	ARG82173 Ketone Body Assay Kit (Colorimetric) is a detection kit for the quantification of Ketone Body (acetoacetate and 3-hydroxybutyrate) in Serum, plasma and urine. This kit is sufficient for up to 100 assays of 3-hydroxybutyric acid and up to 100 assays of acetoacetic acid in 96 well plates.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Ketone Body
Conjugation Note	Read at 340 nm
Sensitivity	0.12 mM
Detection Range	0.12 - 8 mM
Sample Type	Serum, plasma and urine
Sample Volume	10 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	30 min

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82173 Ketone Body Assay Kit (Colorimetric) typical data demonstration image

ARG82173 Ketone Body Assay Kit (Colorimetric) results of a typical data with optical density reading at 340 nm.