

ARG82174
Lactose Assay KitPackage: 100 tests
Store at: -20°C**Summary**


| | |
|---------------------|---|
| Product Description | ARG82174 Lactose Assay Kit is a detection kit for the quantification of Lactose in milk and food. |
| Tested Reactivity | Hu, Ms, Rat, All |
| Tested Application | FuncSt |
| Target Name | Lactose |
| Conjugation Note | Read at 570 nm (colorimetric), or 530/585 nm (fluorometric) |
| Sensitivity | 6 µM |
| Detection Range | 17 - 2000 µM |
| Sample Type | milk and food |
| Sample Volume | 20 µl |

Application Instructions

| | |
|------------------|--|
| Application Note | Please note that this kit does not include a microplate. |
| Assay Time | 30 min |

Properties

| | |
|---------------------|---|
| Form | Liquid |
| Storage instruction | Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |



ARG82174 Lactose Assay Kit typical data demonstration image

ARG82174 Lactose Assay Kit results of a typical data with optical density reading at 570 nm.