

ARG82189
Peroxide Assay Kit (Colorimetric)Package: 250 tests
Store at: 4°C

Summary


Product Description	ARG82189 Peroxide Assay Kit (Colorimetric) is a detection kit for the quantification of Peroxide in serum, citrate-plasma, urine, cell lysate and cell culture supernatants.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Peroxide
Conjugation Note	Read at 585 nm
Sensitivity	0.2 µM
Detection Range	0.2 µM (7 ng/mL) - 30 µM (1,020 ng/mL) H ₂ O ₂
Sample Type	serum, citrate-plasma, urine, cell lysate and cell culture supernatants
Sample Volume	40 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	30 min

Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82189 Peroxide Assay Kit (Colorimetric) typical data demonstration image

ARG82189 Peroxide Assay Kit (Colorimetric) results of a typical data with optical density reading at 585 nm.