

ARG82196 Urease Activity Assay Kit

Package: 100 tests
Store at: 4°C

Summary


Product Description	ARG82196 Urease Activity Assay Kit can be used to measure Urease activity in biological and environmental samples.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Urease
Conjugation Note	Read at 670 nm
Sensitivity	0.003 U/L
Sample Type	biological and environmental samples
Sample Volume	90 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	40 min

Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82196 Urease Activity Assay Kit typical data demonstration image

ARG82196 Urease Activity Assay Kit results of a typical data with optical density reading at 670 nm.