

ARG82218 8-iso-Prostaglandin F2 alpha ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG82218 8-iso-Prostaglandin F2 alpha ELISA Kit is an Enzyme Immunoassay kit for the quantification of 8-iso-Prostaglandin F2 alpha in serum, plasma, urine, cell lysate and tissue lysate.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	8-iso-Prostaglandin F2 alpha
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sample Type	Serum, plasma, urine, cell lysate and tissue lysate.
Standard Range	49 - 200000 pg/ml
Sample Volume	55 µl

Application Instructions


Assay Time	~ 3 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
-----------	--



ARG82218 8-iso-Prostaglandin F2 alpha ELISA Kit standard curve image

ARG82218 8-iso-Prostaglandin F2 alpha ELISA Kit results of a typical standard run with optical density reading at 450 nm.