

ARG82237 MuLV Core Antigen / MuLV p30 ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG82237 MuLV Core Antigen / MuLV p30 ELISA Kit is an Enzyme Immunoassay kit for the quantification of MuLV Core Antigen in cell culture medium.
Tested Reactivity	Virus
Tested Application	ELISA
Target Name	MuLV Core Antigen / MuLV p30
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	300 pg/ml
Sample Type	Cell culture medium
Standard Range	313 - 20000 pg/ml
Sample Volume	100 µl

Application Instructions


Assay Time	4.5 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	New ELISA data calculation tool: Simplify the ELISA analysis by GainData
-----------	---



ARG82237 MuLV Core Antigen / MuLV p30 ELISA Kit standard curve image

ARG82237 MuLV Core Antigen / MuLV p30 ELISA Kit results of a typical standard run with optical density reading at 450 nm.