

ARG82778 Glutamine Synthetase Activity Assay Kit (Colorimetric)

Package: 96 wells
Store at: 4°C and -20°C

Summary

Product Description	ARG82778 Glutamine Synthetase Activity Assay Kit (Colorimetric) is a detection kit for the quantification of Glutamine Synthetase Activity in tissue extract and cell lysate.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Glutamine Synthetase
Conjugation Note	Read at 540 nm.
Detection Range	0.05 - 5 mmol/L
Sample Type	Tissue extract and cell lysate.
Sample Volume	35 µl


Application Instructions

Assay Time	30 min
------------	--------

Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images



ARG82778 Glutamine Synthetase Activity Assay Kit (Colorimetric) standard curve image

ARG82778 Glutamine Synthetase Activity Assay Kit (Colorimetric) results of a typical standard run with optical density reading at 540 nm.