

**ARG82948**  
**Monkey IgA ELISA kit**Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG82948 Monkey IgA ELISA kit is an Enzyme Immunoassay kit for the quantification of Monkey IgA in serum, plasma and cell culture supernatants.
Tested Reactivity	Mk
Tested Application	ELISA
Target Name	IgA
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.78 ng/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	1.56 - 100 ng/ml
Sample Volume	100 µl
Precision	Intra-Assay CV: less than 10% Inter-Assay CV: less than 10%

### Application Instructions

Assay Time	~ 3.5 hours
------------	-------------

### Properties

Form	96 well
Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.