

**ARG83087**  
**IgA ELISA Kit**Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG83087 IgA ELISA Kit is an Enzyme immunoassay kit for the quantitative of IgA ELISA Kit
Tested Reactivity	Other
Tested Application	ELISA
Specificity	Not reacts with IgG and IgM.
Target Name	IgA
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	1.7 µg/ml
Sample Type	Saliva
Standard Range	6.9 – 400 µg/ml
Sample Volume	25 µl
Alternate Names	Immunoglobulin A

### Application Instructions

Assay Time	~2 hours
------------	----------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.