

ARG83088 Fish Histamine ELISA Kit (Fast Track)

Package: 48 wells
Store at: 4°C

Summary

Product Description	ARG83088 Fish Histamine ELISA Kit (Fast Track) is an Enzyme rapid semi-quantitative kit for the quantitative of histamine ELISA Kit
Tested Reactivity	Fsh
Tested Application	ELISA
Target Name	Histamine
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.45 ppm
Standard Range	3.0 – 300 ppm
Sample Volume	50 µl
Alternate Names	Histamine

Application Instructions

Assay Time	~0.5 hour
------------	-----------

Properties

Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.