

ARG83101 Blood / Tissue / Cell Genomic DNA Extraction Kit

Package: 100 tests
Store at: 4°C

Summary

Product Description	Blood / Tissue / Cell Genomic DNA Extraction Kit is designed for the simple and efficient extraction of genome. Blood / Tissue / Cell Genomic DNA Extraction Kit can be applied to extract a small quantity of samples manually and perform in a high-throughput scale automatically. Genomic DNA Extracted by this kit can be used to detect host cell DNA in some experiments.
Tested Application	FuncSt


Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related news: Quality control for engineered T cell therapy
-----------	--

Images



ARG83101 Blood/Tissue/Cell Genomic DNA Extraction Kit in 1% agarose gels

ARG83101 Blood/Tissue/Cell Genomic DNA Extraction Kit with Blood samples.