

ARG83555
Hydroxyl Radical Assay KitPackage: 100 tests
Store at: 4°C

Summary

Product Description	ARG83555 Hydroxyl Radical Assay Kit can be used to measure Hydroxyl Radical in Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Hydroxyl Radical
Conjugation Note	Read at 510 nm
Sensitivity	0.1 µmol/ml
Sample Type	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Standard Range	0.2 µmol/ml - 20 µmol/ml
Sample Volume	10 µl
Alternate Names	Hydroxyl Free Radical

Application Instructions

Assay Time	5 min
------------	-------

Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.