

ARG83557
CA / Carbonic Anhydrase Activity Assay Kit

Package: 100 tests

Store at: 4°C

Summary

Product Description	ARG83557 CA / Carbonic Anhydrase Activity Assay Kit can be used to measure CA / Carbonic Anhydrase in Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	CA / Carbonic Anhydrase
Conjugation Note	Read at 405 nm
Sensitivity	10 mmol/L
Sample Type	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Standard Range	20 mmol/L - 2000 nmol/ml
Sample Volume	20 µl
Alternate Names	EC 4.2.1.1; carbonate dehydratase; carbonate hydro-lyase; carbonic dehydratase

Application Instructions

Assay Time	5 min
------------	-------

Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.