

## ARG83558 GAPDH Activity Assay Kit

Package: 100 tests  
Store at: 4°C

### Summary

Product Description	ARG83558 GAPDH Activity Assay Kit can be used to measure GAPDH in Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	GAPDH
Conjugation Note	Read at 492 nm
Sensitivity	2 nmol/mL
Sample Type	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Standard Range	4 nmol/mL - 400 nmol/mL
Sample Volume	10 µl
Alternate Names	EC 1.2.1.12; glyceraldehyde-3-phosphate dehydrogenase; NAD-dependent glyceraldehyde-3-phosphate dehydrogenase; triosephosphate dehydrogenase

### Application Instructions

Assay Time	5 min
------------	-------

### Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.