

**ARG83559**  
**ADA Activity Assay Kit**Package: 100 tests  
Store at: 4°C

### Summary

Product Description	ARG83559 ADA Activity Assay Kit can be used to measure ADA in Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	ADA
Conjugation Note	Read at 620 nm
Sensitivity	0.005 µmol/ml
Sample Type	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Standard Range	0.01 µmol/ml - 1 µmol/ml
Sample Volume	10 µl
Alternate Names	EC 3.5.4.4; Adenosine aminohydrolase

### Application Instructions

Assay Time	10 min
------------	--------

### Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.