

ARG83561
alpha HBDH Activity Assay KitPackage: 100 tests
Store at: 4°C

Summary

Product Description	ARG83561 alpha HBDH Activity Assay Kit can be used to measure alpha HBDH in Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	alpha HBDH
Conjugation Note	Read at 450 nm
Sensitivity	1.5 nmol/mL
Sample Type	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Standard Range	3 nmol/mL - 300 nmol/mL
Sample Volume	10 µl
Alternate Names	a-HBDH, alpha-HBDH

Application Instructions

Assay Time	5 min
------------	-------

Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.