

ARG83568 AHCY / Adenosylhomocysteinase Activity Assay Kit

Package: 100 tests

Store at: 4°C

Summary

Product Description	ARG83568 AHCY / Adenosylhomocysteinase Activity Assay Kit can be used to measure AHCY / Adenosylhomocysteinase in Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	AHCY / Adenosylhomocysteinase
Conjugation Note	Read at 620 nm
Sensitivity	0.005 µmol/ml
Sample Type	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Standard Range	0.01 - 1 µmol/ml
Sample Volume	10 µl
Alternate Names	EC 3.13.2.1; adenosylhomocysteine hydrolase; AdoHcyase; S-adenosylhomocysteinase; S-adenosylhomocysteine hydrolase; S-adenosylhomocysteine synthase; SAHase

Application Instructions

Assay Time	10 min
------------	--------

Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.