

**ARG83593**  
**Guanosine Assay Kit**Package: 100 assay  
Store at: 4°C, -20°C**Summary**

Product Description	ARG83593 Guanosine Assay Kit is an Enzyme Immunoassay kit for the quantification of Guanosine.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Guanosine
Conjugation	HRP
Conjugation Note	Read at 590 nm.
Sensitivity	0.79 µM
Sample Type	Cell lysate, Tissue lysate and Cell culture supernatant.
Standard Range	1.57 - 100 µM
Sample Volume	50 µL

**Properties**

Form	96 well
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.