

ARG83597
DDDDK tag ELISA KitPackage: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83597 DDDDK tag ELISA Kit is an Enzyme Immunoassay kit for the quantification of DDDDK tag.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	DDDDK tag
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	< 0.05 µg/mL
Sample Type	Cell lysate, Tissue lysate and Purified protein
Standard Range	0.078 - 5 µg/mL
Sample Volume	50 µL

Properties

Form	96 well
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.