

ARG83598
E. coli O157 ELISA KitPackage: 96 wells
Store at: 4°C, -80°C

Summary

Product Description	ARG83598 E. coli O157 ELISA Kit is an Enzyme Immunoassay kit for the quantification of E. coli O157.
Tested Reactivity	E. coli
Tested Application	ELISA
Target Name	E. coli O157
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sample Type	Food, Fecal and Environmental Sample
Standard Range	10^4 CFU equivalents/mL - 5×10^7 CFU equivalents /mL
Sample Volume	100 μ L

Properties

Form	96 well
Storage instruction	Store the kit at 4°C, -80°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.