

**ARG83599**  
**Anoikis Assay Kit**Package: 96 wells  
Store at: 4°C, -20°C**Summary**

Product Description	ARG83599 Anoikis Assay Kit is a detection kit for the quantification of Anoikis in cell.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Anoikis
Conjugation Note	Read at 485/525 nm.
Sample Type	Cell culture supernatant.

**Properties**

Form	96 well
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.