

ARG83600 Cytotoxicity Assay Kit

Package: 96 assay
Store at: 4°C, -20°C

Summary

Product Description	ARG83600 Cytotoxicity Assay Kit is a detection kit for the quantification and monitoring of cell viability and growth.
Tested Reactivity	Other
Tested Application	FuncSt
Specificity	The kit does not react with bacteria or yeast.
Target Name	Cytotoxicity
Conjugation Note	Read at 570 nm.
Sample Type	Cell culture supernatant.

Properties

Form	96 well
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.