

ARG83605
Adenosine Monophosphate Assay KitPackage: 100 assay
Store at: 4°C, -20°C**Summary**

Product Description	ARG83605 Adenosine Monophosphate Assay Kit is an detection kit for the quantification of Adenosine Monophosphate.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Adenosine
Conjugation Note	Read at 540 nm.
Sensitivity	0.78 μ M
Sample Type	Serum, plasma, urine, saliva, cell/tissue lysate and cell culture supernatants.
Standard Range	1.56 - 100 μ M
Sample Volume	50 μ L

Properties

Form	96 well
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.