

ARG83606 Phosphoenolpyruvate Assay Kit

Package: 100 assay
Store at: 4°C, -20°C

Summary

Product Description	ARG83606 Phosphoenolpyruvate Assay Kit is a detection kit for the quantification of Phosphoenolpyruvate.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Phosphoenolpyruvate
Conjugation Note	Read at 540 nm.
Sensitivity	1.57 μ M
Sample Type	Serum, plasma, urine and cell culture supernatants.
Standard Range	3.13 - 200 μ M
Sample Volume	50 μ L

Properties

Form	96 well
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.