

ARG83608
Tryptophan Assay KitPackage: 100 assay
Store at: 4°C, -20°C**Summary**

Product Description	ARG83608 Tryptophan Assay Kit is a detection kit for the quantification of Tryptophan.
Tested Application	FuncSt
Target Name	Tryptophan
Conjugation Note	Read at 540 nm.
Sensitivity	0.78 µM
Sample Type	Serum, plasma, urine and cell culture supernatants.
Standard Range	1.56 - 100 µM
Sample Volume	50 µL

Properties

Form	96 well
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.