

**ARG83613**  
**ROCK Activity ELISA Kit**Package: 96 wells  
Store at: 4°C, -20°C

### Summary

Product Description	ROCK Activity ELISA Kit is an Enzyme Immunoassay kit for the quantification of ROCK Activity in purified kinase, cell/tissue lysate.
Tested Reactivity	All
Tested Application	ELISA
Target Name	ROCK Activity
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sample Type	Tissue/cell lysate samples.
Sample Volume	90 µl
Alternate Names	Rho Kinase

### Properties

Form	Liquid
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.