

ARG83617
Zika virus NS1 Protein ELISA KitPackage: 96 wells
Store at: 4°C, -20°C

Summary


Product Description	ARG83617 Zika virus NS1 Protein ELISA Kit is an Enzyme Immunoassay kit for the quantification of Zika virus NS1 in Purified virus or unpurified viral supernatant.
Tested Reactivity	Virus
Tested Application	ELISA
Target Name	Zika virus NS1
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	10 pg/ml
Sample Type	Purified virus or unpurified viral supernatant.
Standard Range	19.5 - 1250 pg/ml
Sample Volume	100 µl

Application Instructions

Assay Time	3.5 hours
------------	-----------

Properties

Form	Liquid
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG83617 Zika virus NS1 Protein ELISA Kit standard curve image

ARG83617 Zika virus NS1 Protein ELISA Kit results of a typical standard run with optical density reading at 450 nm.