

## ARG83622 Influenza B nucleoprotein ELISA Kit

Package: 96 wells  
Store at: 4°C, -20°C

### Summary

Product Description	ARG83622 Influenza B nucleoprotein ELISA Kit is an Enzyme Immunoassay kit for the quantification of Influenza B nucleoprotein in Purified virus or unpurified viral supernatant.
Tested Reactivity	Influenza B virus
Tested Application	ELISA
Target Name	Influenza B nucleoprotein
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	0.15 ng/ml
Sample Type	Purified virus or unpurified viral supernatant.
Standard Range	0.313 - 20 ng/ml
Sample Volume	100 µl

### Application Instructions

Assay Time	4.5 hours
------------	-----------

### Properties

Form	Liquid
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.