

ARG83734
FSME (TBEV) IgM ELISA KitPackage: 96 wells
Store at: 4°C

Summary

Product Description	ARG83734 FSME (TBEV) IgM ELISA Kit is an Enzyme Immunoassay kit for the qualitative of FSME (TBEV) IgM in Serum
Tested Reactivity	Animals
Tested Application	ELISA
Target Name	FSME (TBEV) IgM
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sample Type	Serum
Standard Range	Cut-off
Sample Volume	100 µl

Application Instructions

Assay Time	2 hours
------------	---------

Properties

Form	96 well
Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.